

Cuattro Medical Releases Cuattro **one** Long Format Digital Detector

NEWS PROVIDED BY
Cuattro Medical LLC
July 01, 2019, 09:30 ET

GOLDEN, Colo., July 01, 2019 /PRNewswire/ -- Cuattro Medical LLC, an industry leading developer and provider of innovative digital radiography solutions is proud to introduce the new Cuattro ONE Long Format Detector. The Cuattro ONE is a 17" x 51" long format detector designed to simplify the challenges associated with full spine and long leg imaging studies. A single X-ray exposure produces a full field image preview (17" x 50.6") in 6 seconds with a final, fully processed image in 12 seconds. Additionally, the user can isolate the top, middle or bottom portions of the 51" array to complete conventional 17" x 17" weight bearing views. Available in either floor mounted or mobile wall stand configurations, the Cuattro ONE can be retrofitted to existing X-ray equipment and suites or included as a vital part of a new full room installation. The design and manufacturer of the new ONE was a joint development collaboration by Vieworks Corporation (<https://www.vieworks.com/>) and Cuattro Medical LLC (<https://cuattro.com/>).

For additional details on the new Cuattro ONE, please visit <https://cuattro.com/>

About Cuattro Medical:

Cuattro is a world leader in Digital Radiography system design, installation and support. With over 20 years of experience, the Cuattro leadership team has successfully partnered with doctors and institutions across numerous specialties resulting in over 6,000 successful Flat Panel Digital Radiography installations worldwide. Cuattro's innovative Digital Radiography platforms include both full room and retrofit solutions for Orthopedic, Hospital, Urgent Care and Family Practice applications. Cuattro products are also available for OEM and NDT use.

Contact: Clint Roth, 760 518-5220 croth@cuattro.com