






Experience the Cuattro ONE difference! Apply to qualify for our free, no obligation in-clinic 30 day trial.*

Long format X-ray exams performed with CR and conventional DR with Stitching are difficult and time consuming. The usual result is underutilization of Long Format Studies and frustration with the wait times and results of many of the exams that are completed.



one™

Long Format Digital Radiography

- ONE exposure
ONE image
ONE minute

17x51\"/>



CR Long Format Exam: 15 - 25 Minutes

Computed Radiography long format studies are time consuming to complete, making patients and physicians wait far too long for study results.

Conventional DR Long Format: 8 - 12 Minutes

Conventional DR Long Format Exams require up to three separate exposures, increasing patient and staff x-ray dose exposure. Multiple exposures often result in patient motion, necessitating retakes.

Cuattro ONE Long Format: Under 1 Minute

Cuattro ONE allows the user to complete full length (50.5\") Long Format Exams in less than a minute with a single exposure. This results in lower total x-ray dose exposure to patients and staff.

ONE exposure.
ONE image.
ONE minute.

The speed and efficiency of the Cuattro ONE Long Format Digital Radiography system increases study utilization over 200% on average. Since no stitching is required, clinical accuracy improves with significantly reduced patient motion and much better anatomical alignment. Due to a substantially improved workflow, diagnostic suites have more time to perform more X-ray studies.

www.cuattro.com | toll free: 800.709.4514

*Free 30 day trial offer subject to change or cancellation at any time without notice. Restrictions apply. © 2019 Cuattro, LLC. All rights reserved. Specifications and option availability are subject to change without notice. Consult a Cuattro sales specialist for configuration details. Names of companies or products appearing in this document are trademarks of their respective owners. Cuattro and Cuattro ONE are trademarks of Cuattro, LLC in the United States and may also be registered trademarks or trademarks in other countries. 19LT1001

